

**PRODUCT DATA SHEET**
**TMTsixplex™ Label Reagent Set, 6 Samples**
**Product Number:** 90061

**Lot Number:** TC255606

**Form:** The TMTsixplex™ Label Reagents are supplied dried, 0.8 mg/tube. Make a stock solution by reconstituting each tube with 41 µl acetonitrile.

**Note:** Please refer to Table 2 of the instruction booklet for TMTsixplex mass modification, reporter ion masses or mass tolerance window in your data analysis software.

\*\* Reporter ion isotopic distributions (-2, -1, +1, +2) are primarily for carbon isotopes with reporter ion interference for each mass tag shown. Reporter ion isotopic distributions can be used as isotope correction factors in TMT6plex method template in Proteome Discoverer software version 1.4 and above.

**\*\*Reporter Ion Isotopic Distributions:**

Mass Tag	Mass Reporter	Reporter Ion	-2	-1	Monoisotopic	+1	+2
TMT <sup>6</sup> -126	126	126.127726	0.0%	0.0%	100%	8.2%	0.7%
TMT <sup>6</sup> -127	127N	127.124761	0.0%	0.5%	100%	5.9%	0.0%
TMT <sup>6</sup> -128	128C	128.134436	0.0%	0.9%	100%	6.8%	0.0%
TMT <sup>6</sup> -129	129N	129.131471	0.0%	4.5%	100%	4.8%	0.0%
TMT <sup>6</sup> -130	130C	130.141145	0.5%	3.0%	100%	3.1%	0.0%
TMT <sup>6</sup> -131	131	131.138180	0.0%	3.5%	100%	3.1%	0.0%

**Stability:** One year from date of product receipt.

**Storage:** Store at -20°C.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).** There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

 Pierce Biotechnology  
 3747 N. Meridian Road

 PO Box 117  
 Rockford, IL 61105 USA

 (800) 874-3723 tel  
 (815) 968-0747 tel

 (815) 968-7316 fax  
[www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce)