

PRODUCT DATA SHEET

TMTsixplex[™] Label Reagent Set, 6 Samples

Product Number: 90061 Lot Number: QC213931

Form: The TMTsixplex™ Label Reagents are supplied dried, 2 tubes x 0.8 mg/tube. Make a stock solution by reconstituting each tube with 41 µl acetonitrile.

Note: The positions of ¹⁵N and ¹³C heavy isotopes have changed recently in the TMTsixplex-127 and TMTsixplex-129 chemical structures. TMTsixplex reagents now permit full sixplex quantitation with both CID and ETD fragmentation methods. The 127 and 129 Da reporter ion masses have changed slightly; please refer to Table 3 of the instruction booklet to update the reporter ion masses or mass tolerance window in your data analysis software.

Reporter Ion Isotopic Distributions:

Label Reagent	-2	-1	Main (Normalized)	+1	+2
TMT ⁶ -126	0.0%	0.0%	100%	7.1%	0.0%
TMT ⁶ -127	0.0%	0.4%	100%	5.4%	0.0%
TMT ⁶ -128	0.0%	1.1%	100%	4.9%	0.0%
TMT ⁶ -129	0.0%	1.5%	100%	4.4%	0.0%
TMT ⁶ -130	0.1%	2.9%	100%	2.4%	0.0%
TMT ⁶ -131	0.2%	2.9%	100%	2.5%	0.0%

Stability: One year from date of product receipt.

Storage: Store at -20°C.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").